



SigmaStat (Windows)

Ver.1.0. 89


Installation Guide



SigmaStat 1.0 Installation Guide:

Download Instructions:

1. Uninstall SigmaStat (If previous version is installed – Steps given below).
2. **Remove the folder “/My Documents/SigmaStat/SSNG1/”** (If any previous version is installed), this folder will be created during SigmaStat application startup next time.
3. Download (SigmaStat Setup 1.0.89.zip) file from the given location given.
4. Unzip the SigmaStat windows installer package (SigmaStat Setup 1.0.89.zip) into a separate local folder (Ex. D:\SigmaStat_Installer). Ensure that the **"Libraries"** folder is present in the same directory as the installer during installation as shown below.



5. Install the new build using the installer file **"SigmaStat Setup 1.0.89.exe"** extracted from the ZIP file.


Notes: Please note that by default SigmaStat is installed with 15 days trial license.

I) Installation Steps SigmaStat 1.0 in Windows:


Step 1: Navigate to the folder where “SigmaStat Setup 1.0.89.exe” installer has been extracted.

Step 2: Execute the SigmaPlot ng installer file “SigmaStat Setup 1.0.89.exe”


Step 3: Accept the license agreement as shown below.




Step 4: Choose Install Location




Step 5: Installation starts with progress dialog




Step 6: SigmaStat Ollama Installation will start automatically, if it is not installed before.



Step 7: Finish installation and Run SigmaStat.




Step 8: SigmaStat application starts with default desktop window.




II) License Activation:

Step 1: Run SigmaStat, select thru ribbon “Assist->License Status” to display ‘License Utility’




Step 2: License Utility dialog with annual license




Step 3: Select the license product key or license file:

Select either product key or license file as shown below and apply the license. Product key is supplied as text file and need to copy the product key code in the file. License file is supplied as a file with extension V2C.


Product Key:



License File:





Step 4: Restart the application for the new license:



III) Uninstall SigmaStat 1.0 in Windows:

Use “Add or Remove Programs” to uninstall the following applications.

Uninstall SigmaStat:



Uninstall SigmaStat Ollama:

Apps > Installed apps

 Ollama version 0.31.1
0.31.1 | Ollama | 11-08-2026

3.03 GB ...

Ollama Uninstall

Uninstall Status

Please wait while Ollama is removed from your computer.



Uninstalling Ollama...

Cancel